

## ABC IHC Protocol

8/1/2008

The following outlines a DAB/peroxidase protocol, using avidin-biotin as an amplifying step, on paraffin-embedded tissues, to allow optimal colorimetric detection of markers for immunohistochemistry.

### Materials needed:

- Antibody Amplifier*<sup>TM</sup> (one Antibody Amplifier fits a maximum of 12 slides)
- Amplifying Antibody Dilution Buffer (1x) cat# AA3
- Primary antibody (try to use the same LOT of antibody for all slides in a project)
- Amplifying IHC Wash Buffer (20x) cat # AA4 (you will need to dilute 20 fold in ddH<sub>2</sub>O for working dilution)
- Epitope Unmasking Solution (10x) cat# AA5 (you will need to dilute 10 fold in ddH<sub>2</sub>O for working dilution)
- Vectastain Kit (Vector Labs; contains biotinylated 2°, reagents A and B)
- DAB reagents
- Xylene
- Ethanol
- 1x PBS
- Peroxidase Blocking Solution (3% H<sub>2</sub>O<sub>2</sub> in PBS):
  - 30% H<sub>2</sub>O<sub>2</sub> ----- 2 ml
  - 1X PBS ----- 18 ml
- Coverslips
- Methyl Green (1% in water) or Hematoxylin (Gill-1)
- Slide racks and holder:
- Permount
- ddH<sub>2</sub>O

**Approximate time involved: 4 hours day 1; 4-5 hours day 2;**

### Procedure:

1. Put slides in a glass slide holder and place in the following solutions serially (at RT):
  - a. xylene: 5 min., swishing every 2 min. (IN FUME HOOD)
  - b. xylene: 5min., swishing every 2 min. (IN FUME HOOD)
  - c. 100% ETOH: 2 min., swishing every 30 sec.
  - d. 95% ETOH: 2 min., rocking.
  - e. 80% ETOH: 2 min., rocking.
  - f. 70% ETOH: 2min., rocking.
  - g. 50% ETOH: 2 min., rocking.
  - h. dd H<sub>2</sub>O: 2 X 3min., rocking.
2. Wash with 1X Amplifying IHC Wash Buffer (cat # AA4) 2 x 5 min., rocking at RT.
3. Wash 2 X 3 min. in 1X Amplifying IHC Wash Buffer (cat # AA4)
4. Transfer slide holder to plastic, autoclavable container and immerse in 1X Epitope Unmasking Solution (cat# AA5).

5. Autoclave for 20-30 min., while covered with foil.

6. Retrieve sample from autoclave and place on counter to let cool at RT 20 min. (continuing emersion in Epitope Unmasking Solution (cat #AA5).

NOTE: for fatty tissue, such as brain or breast, use these directions for antigen retrieval: place slides in slide rack. Lower slide rack into gently boiling 1X Epitope Unmasking Solution (cat# AA5) (about 600 mls of buffer). We use a large beaker on a heat block. Boil gently for exactly 8 min. Remove from heat source. Slowly add 500 mls dd H<sub>2</sub>O to bring buffer temp down. Quickly put slide rack in glass slide holder and wash 2 X 5 min in 1X Amplifying IHC Wash Buffer (cat # AA4)

7. Wash with 1X Amplifying IHC Wash Buffer (cat # AA4) 2 x 5 min., rocking at RT.

8. Wipe off back and sides of slide with a kimwipe, lay flat on benchtop (in a humidity chamber if desired) and cover tissue with 1-3 drops of peroxidase block (3% H<sub>2</sub>O<sub>2</sub> in PBS). Incubate RT 30 min., checking every 5-10 min. to make sure samples don't dry out. If starting to dry, add more peroxidase block.

9. Put slides back into slide carrier and wash 3 X 3min. in 1X Amplifying IHC Wash Buffer (cat # AA4)

10. Dilute 1° antibody in Amplifying Antibody Dilution Buffer (cat# AA3). Use a minimum of 3 mls per slide for Antibody Amplifier. The dilution will, of course, vary for each antibody, but in our experience, you can use at least 5x and often >20x the dilution you typically use for the standard sedentary method. Place lid onto the *Antibody Amplifier<sup>TM</sup>*, and place on a standard laboratory rocker. Rock overnight at 4°C.

After use, store primary antibody at -20°C.

12. Put slides back into glass slide carrier and wash 2 X 5 min. in 1X Amplifying IHC Wash Buffer (cat # AA4).

13. Incubate slides in biotinylated secondary antibody. Secondary is provided as part of the Vectastain kit by Vector Labs. Secondary antibody dilution is 1:2000 in Amplifying Antibody Dilution Buffer (cat# AA3). Incubate slides rocking at RT 1 hour in *Antibody Amplifier<sup>TM</sup>*. Remember that each slide needs 3 mls.

14. As soon as slides are incubating in secondary antibody, mix up an ABC solution (found in Vectastain kit). The ABC solution is diluted 1:1000. Remember if using *Antibody Amplifier<sup>TM</sup>*, to make 3mls per slide. Let the mixed solution sit RT 30-60 min. before using.

15. Put slides in glass slide carrier and wash 2 X 5 min. in 1X Amplifying IHC Wash Buffer (cat # AA4).

**16.** Incubate slides in premade ABC solution, rocking at RT 1 hr in the *Antibody Amplifier<sup>TM</sup>*.

**17.** Make up DAB solution according to manufacturer's instructions. You need 100-200 ul per slide. Incubate RT 20-30 min.

**18.** Wash 2 X 5min in 1X Amplifying IHC Wash Buffer (cat # AA4).

**19.** Wash 1 X 3 min dd H<sub>2</sub>O.

**20.** Wipe slides around tissue with kimwipe, keeping tissue moist, and add DAB solution. Check the staining of one slide first by testing a positive control: Incubate with 100-200 ul of DAB solution, starting at 30 seconds. Check using the microscope at 40X. Keep note of DAB times. Wash slides as soon as staining is optimal with ddH<sub>2</sub>O.

**NOTE:** All slides that you are comparing must be incubated for the same time period with DAB once this time is established with the positive control.

**21.** Wash slides in ddH<sub>2</sub>O 1 X 3 min.

**22.** Counterstain using Methyl Green: Filter 3-10mls of methyl green using a syringe filter into a 15 ml tube. Add 200 ul of methyl green onto tissue and incubate at RT 2 min (or whatever time you have established). Rinse in containers of ddH<sub>2</sub>O until slides are clean and let air-dry, with slides standing vertically on their ends.

**23.** After tissues are completely dry, add one drop of Permount, put on coverslip, and smooth out bubbles using big end of pipettor tip. Label slides with: date, antibody, dilution, and timing. After 4-5 hours, clean off extra permount with xylene (xylene is optional-only if covering tissue).